

Compendium of overhangs

There are **136** overhangs in this report.

Please see the Appendix on the last page for an explanation of details.



GTTT

GC content: 25 %.

Has 3 identical bases in a row. However, this has not shown to be very important.

Can form the following amino acids in 6 translation frames:

K[PRQLH]
[PEGRSA*IKTQVL]N
[EQK*][T]
V[FWSCY*L]
[PEGWRSA*KTQMVL]F
[SGCR][LF]



TGTT

GC content: 25 %.

Can form the following amino acids in 6 translation frames:

N[RSIKTNM]

[PEGRSA*IKTQVL]T

[EQK*][HQ]

C[FWSCY*L]

[PGFSRACYIDTNVLH]V

[MVL][LF]



GGTT

GC content: 50 %.

Can form the following amino acids in 6 translation frames:

N[PRQLH]
[PEGRSA*IKTQVL]T
[EQK*][P]
G[FWSCY*L]
[PEGWRSA*KTQMVL]V
[RGW][LF]

Page 4 of 138



AGTT

GC content: 25 %.

Can form the following amino acids in 6 translation frames:

N[FWSCY*L]
[PEGRSA*IKTQVL]T
[EQK*][L]
S[FWSCY*L]
[PEGRSA*IKTQVL]V
[EQK*][LF]



ACTT

GC content: 25 %.

Can form the following amino acids in 6 translation frames:

K[FWSCY*L]
[PEGRSA*IKTQVL]S
[EQK*][V]
T[FWSCY*L]
[PEGRSA*IKTQVL]L

[HDYN][LF]



TATT

GC content: 0 %.

Can form the following amino acids in 6 translation frames:

N[RSIKTNM]
[PEGRSA*IKTQVL]I
[EQK*][Y*]
Y[FWSCY*L]
[PGFSRACYIDTNVLH]I

[IVL][LF]



GATT

GC content: 25 %.

Can form the following amino acids in 6 translation frames:

N[PRQLH]
[PEGRSA*IKTQVL]I
[EQK*][S]
D[FWSCY*L]
[PEGWRSA*KTQMVL]I
[RG*][LF]



AATT

GC content: 0 %.

The overhang is palindromic, cannot be used for DNA assembly.

Can form the following amino acids in 6 translation frames:

N[FWSCY*L]

[PEGRSA*IKTQVL]I

[EQK*][LF]

N[FWSCY*L]

[PEGRSA*IKTQVL]I

[EQK*][LF]



TTGT

GC content: 25 %.

Can form the following amino acids in 6 translation frames:

T[RSIKTNM]
[PEGRSA*IKTQVL]Q
[HDYN][KN]
L[FWSCY*L]
[PGFSRACYIDTNVLH]C

[IVFL][V]



CTGT

GC content: 50 %.

Can form the following amino acids in 6 translation frames:

T[GEADV]
[PEGRSA*IKTQVL]Q
[HDYN][RS]
L[FWSCY*L]
[PGFSRACYIDTNVLH]C
[SPAT][V]



TGGT

GC content: 50 %.

Can form the following amino acids in 6 translation frames:

T[RSIKTNM]
[PEGRSA*IKTQVL]P
[HDYN][HQ]
W[FWSCY*L]
[PGFSRACYIDTNVLH]G
[MVL][V]

Page 12 of 138



CGGT

GC content: 75 %.

[SPAT][V]

Can form the following amino acids in 6 translation frames:

T[GEADV]
[PEGRSA*IKTQVL]P
[HDYN][R]
R[FWSCY*L]
[PGFSRACYIDTNVLH]G



TCGT

GC content: 50 %.

[IVFL][V]

Can form the following amino acids in 6 translation frames:

T[RSIKTNM]
[PEGRSA*IKTQVL]R
[HDYN][ED]
S[FWSCY*L]
[PGFSRACYIDTNVLH]R

Page 14 of 138



GCGT

GC content: 75 %.

Can form the following amino acids in 6 translation frames:

T[PRQLH]
[PEGRSA*IKTQVL]R
[HDYN][A]
A[FWSCY*L]
[PEGWRSA*KTQMVL]R

[SGCR][V]



ACGT

GC content: 50 %.

The overhang is palindromic, cannot be used for DNA assembly.

Can form the following amino acids in 6 translation frames:

T[FWSCY*L]
[PEGRSA*IKTQVL]R
[HDYN][V]
T[FWSCY*L]
[PEGRSA*IKTQVL]R
[HDYN][V]



TAGT

GC content: 25 %.

[IVL][V]

Can form the following amino acids in 6 translation frames:

T[RSIKTNM]
[PEGRSA*IKTQVL]L
[HDYN][Y*]
*[FWSCY*L]
[PGFSRACYIDTNVLH]S



CAGT

GC content: 50 %.

Can form the following amino acids in 6 translation frames:

T[GEADV]
[PEGRSA*IKTQVL]L
[HDYN][WC*]
Q[FWSCY*L]
[PGFSRACYIDTNVLH]S
[SPAT][V]



GGCT

GC content: 75 %.

Can form the following amino acids in 6 translation frames:

S[PRQLH]
[PEGRSA*IKTQVL]A
[EQK*][P]
G[FWSCY*L]
[PEGWRSA*KTQMVL]A
[RGW][L]



CGCT

GC content: 75 %.

Can form the following amino acids in 6 translation frames:

S[GEADV]

[PEGRSA*IKTQVL]A

[EQK*][R]

R[FWSCY*L]

[PGFSRACYIDTNVLH]A

[SPAT][L]



X AGCT

AGCT

GC content: 50 %.

The overhang is palindromic, cannot be used for DNA assembly.

Can form the following amino acids in 6 translation frames:

S[FWSCY*L]

[PEGRSA*IKTQVL]A

[EQK*][L]

S[FWSCY*L]

[PEGRSA*IKTQVL]A

[EQK*][L]



TCCT

GC content: 50 %.

[IVFL][L]

Can form the following amino acids in 6 translation frames:

R[RSIKTNM]
[PEGRSA*IKTQVL]G
[EQK*][ED]
S[FWSCY*L]
[PGFSRACYIDTNVLH]P





GC content: 75 %.

Has 3 identical bases in a row. However, this has not shown to be very important.

Can form the following amino acids in 6 translation frames:

R[GEADV]
[PEGRSA*IKTQVL]G
[EQK*][G]
P[FWSCY*L]
[PGFSRACYIDTNVLH]P
[SPAT][L]



AGGT

GC content: 50 %.

Can form the following amino acids in 6 translation frames:

T[FWSCY*L]

[PEGRSA*IKTQVL]P

[HDYN][L]

R[FWSCY*L]

[PEGRSA*IKTQVL]G

[EQK*][V]



TACT

GC content: 25 %.

[IVL][L]

Can form the following amino acids in 6 translation frames:

S[RSIKTNM]
[PEGRSA*IKTQVL]V
[EQK*][Y*]
Y[FWSCY*L]
[PGFSRACYIDTNVLH]T



CACT

GC content: 50 %.

Can form the following amino acids in 6 translation frames:

S[GEADV]

[PEGRSA*IKTQVL]V

 $[EQK^*][WC^*]$

H[FWSCY*L]

[PGFSRACYIDTNVLH]T

[SPAT][L]



ATAT

GC content: 0 %.

The overhang is palindromic, cannot be used for DNA assembly.

Can form the following amino acids in 6 translation frames:

I[FWSCY*L]
[PEGRSA*IKTQVL]Y
[HDYN][MI]
I[FWSCY*L]
[PEGRSA*IKTQVL]Y

[HDYN][MI]



GGAT

GC content: 50 %.

Can form the following amino acids in 6 translation frames:

I[PRQLH]
[PEGRSA*IKTQVL]S
[HDYN][P]
G[FWSCY*L]
[PEGWRSA*KTQMVL]D
[RGW][MI]



CGAT

GC content: 50 %.

Can form the following amino acids in 6 translation frames:

I[GEADV]
[PEGRSA*IKTQVL]S
[HDYN][R]
R[FWSCY*L]
[PGFSRACYIDTNVLH]D

[SPAT][MI]



ATCT

GC content: 25 %.

Can form the following amino acids in 6 translation frames:

R[FWSCY*L]

[PEGRSA*IKTQVL]D

[EQK*][MI]

I[FWSCY*L]

[PEGRSA*IKTQVL]S

[HDYN][L]



TCAT

GC content: 25 %.

The overhang contains the start codon ATG.

The overhang contains a stop codon (TAA, TAG or TGA).

Can form the following amino acids in 6 translation frames:

M[RSIKTNM]
[PEGRSA*IKTQVL]*
[HDYN][ED]
S[FWSCY*L]
[PGFSRACYIDTNVLH]H
[IVFL][MI]



ATGT

GC content: 25 %.

The overhang contains the start codon ATG.

Can form the following amino acids in 6 translation frames:

T[FWSCY*L]
[PEGRSA*IKTQVL]H
[HDYN][MI]
M[FWSCY*L]
[PEGRSA*IKTQVL]C
[HDYN][V]



TAAT

GC content: 0 %.

[IVL][MI]

Can form the following amino acids in 6 translation frames:

I[RSIKTNM]
[PEGRSA*IKTQVL]L
[HDYN][Y*]
*[FWSCY*L]
[PGFSRACYIDTNVLH]N

Page 33 of 138



GAAT

GC content: 25 %.

Can form the following amino acids in 6 translation frames:

I[PRQLH] [PEGRSA*IKTQVL]F

[HDYN][S] E[FWSCY*L]

[PEGWRSA*KTQMVL]N

[RG*][MI]





GC content: 0 %.

Has 3 identical bases in a row. However, this has not shown to be very important.

Can form the following amino acids in 6 translation frames:

K[FWSCY*L]
[PEGRSA*IKTQVL]N
[EQK*][MI]
I[FWSCY*L]
[PEGRSA*IKTQVL]F
[HDYN][LF]



TTTG

GC content: 25 %.

Has 3 identical bases in a row. However, this has not shown to be very important.

Can form the following amino acids in 6 translation frames:

Q[RSIKTNM]
[PGFSRACYIDTNVLH]K
[SPAT][KN]
F[GEADV]
[PGFSRACYIDTNVLH]L
[IVFL][WC*]



GTTG

GC content: 50 %.

Can form the following amino acids in 6 translation frames:

Q[PRQLH]
[PGFSRACYIDTNVLH]N
[SPAT][T]
V[GEADV]
[PEGWRSA*KTQMVL]L
[SGCR][WC*]



CAAT

GC content: 25 %.

Can form the following amino acids in 6 translation frames:

I[GEADV]
[PEGRSA*IKTQVL]L
[HDYN][WC*]
Q[FWSCY*L]
[PGFSRACYIDTNVLH]N
[SPAT][MI]



TGTG

GC content: 50 %.

Can form the following amino acids in 6 translation frames:

H[RSIKTNM]
[PGFSRACYIDTNVLH]T
[SPAT][HQ]
C[GEADV]
[PGFSRACYIDTNVLH]V

[MVL][WC*]



GGTG

GC content: 75 %.

[RGW][WC*]

Can form the following amino acids in 6 translation frames:

H[PRQLH]
[PGFSRACYIDTNVLH]T
[SPAT][P]
G[GEADV]
[PEGWRSA*KTQMVL]V



CGTG

GC content: 75 %.

Can form the following amino acids in 6 translation frames:

H[GEADV]

[PGFSRACYIDTNVLH]T

[SPAT][R]

R[GEADV]

[PGFSRACYIDTNVLH]V

[SPAT][WC*]



TCTG

GC content: 50 %.

Can form the following amino acids in 6 translation frames:

Q[RSIKTNM]

[PGFSRACYIDTNVLH]R

[SPAT][ED]

S[GEADV]

[PGFSRACYIDTNVLH]L

[IVFL][WC*]



GATG

GC content: 50 %.

[RG*][WC*]

The overhang contains the start codon ATG.

Can form the following amino acids in 6 translation frames:

H[PRQLH]
[PGFSRACYIDTNVLH]I
[SPAT][S]
D[GEADV]
[PEGWRSA*KTQMVL]M



CATG

GC content: 50 %.

The overhang is palindromic, cannot be used for DNA assembly.

The overhang contains the start codon ATG.

Can form the following amino acids in 6 translation frames:

H[GEADV]
[PGFSRACYIDTNVLH]M
[SPAT][WC*]
H[GEADV]
[PGFSRACYIDTNVLH]M
[SPAT][WC*]



CATT

GC content: 25 %.

The overhang contains the start codon ATG.

Can form the following amino acids in 6 translation frames:

N[GEADV]
[PEGRSA*IKTQVL]M
[EQK*][WC*]
H[FWSCY*L]
[PGFSRACYIDTNVLH]I

[SPAT][LF]



GTGG

GC content: 75 %.

Can form the following amino acids in 6 translation frames:

P[PRQLH]

[PGFSRACYIDTNVLH]H

[SPAT][T]

V[GEADV]

[PEGWRSA*KTQMVL]W

[SGCR][G]



CCAT

GC content: 50 %.

[SPAT][MI]

The overhang contains the start codon ATG.

Can form the following amino acids in 6 translation frames:

M[GEADV]
[PEGRSA*IKTQVL]W
[HDYN][G]
P[FWSCY*L]
[PGFSRACYIDTNVLH]H



TGGG

GC content: 75 %.

Has 3 identical bases in a row. However, this has not shown to be very important.

Can form the following amino acids in 6 translation frames:

P[RSIKTNM]
[PGFSRACYIDTNVLH]P
[SPAT][HQ]
W[GEADV]
[PGFSRACYIDTNVLH]G
[MVL][G]



GCGG

GC content: 100 %.

Can form the following amino acids in 6 translation frames:

P[PRQLH]

[PGFSRACYIDTNVLH]R

[SPAT][A]

A[GEADV]

[PEGWRSA*KTQMVL]R

[SGCR][G]



X CCGG

CCGG

GC content: 100 %.

The overhang is palindromic, cannot be used for DNA assembly.

Can form the following amino acids in 6 translation frames:

P[GEADV]

[PGFSRACYIDTNVLH]R

[SPAT][G]

P[GEADV]

[PGFSRACYIDTNVLH]R

[SPAT][G]



CCGT

GC content: 75 %.

[SPAT][V]

Can form the following amino acids in 6 translation frames:

T[GEADV]
[PEGRSA*IKTQVL]R
[HDYN][G]
P[FWSCY*L]
[PGFSRACYIDTNVLH]R



GAGG

GC content: 75 %.

Can form the following amino acids in 6 translation frames:

P[PRQLH]

[PGFSRACYIDTNVLH]L

[SPAT][S]

E[GEADV]

[PEGWRSA*KTQMVL]R

[RG*][G]



CCTG

GC content: 75 %.

Can form the following amino acids in 6 translation frames:

Q[GEADV]

[PGFSRACYIDTNVLH]R

[SPAT][G]

P[GEADV]

[PGFSRACYIDTNVLH]L

[SPAT][WC*]



CCTT

GC content: 50 %.

Can form the following amino acids in 6 translation frames:

K[GEADV]

[PEGRSA*IKTQVL]R

[EQK*][G]

P[FWSCY*L]

[PGFSRACYIDTNVLH]L

[SPAT][LF]



GTCG

GC content: 75 %.

Can form the following amino acids in 6 translation frames:

R[PRQLH]

[PGFSRACYIDTNVLH]D

[SPAT][T]

V[GEADV]

[PEGWRSA*KTQMVL]S

[SGCR][R]



TGCG

GC content: 75 %.

Can form the following amino acids in 6 translation frames:

R[RSIKTNM]

[PGFSRACYIDTNVLH]A

[SPAT][HQ]

C[GEADV]

[PGFSRACYIDTNVLH]A

[MVL][R]



GGCG

GC content: 100 %.

Can form the following amino acids in 6 translation frames:

R[PRQLH]
[PGFSRACYIDTNVLH]A
[SPAT][P]
G[GEADV]

[PEGWRSA*KTQMVL]A

[RGW][R]



X CGCG

CGCG

GC content: 100 %.

The overhang is palindromic, cannot be used for DNA assembly.

Can form the following amino acids in 6 translation frames:

R[GEADV]

[PGFSRACYIDTNVLH]A

[SPAT][R]

R[GEADV]

[PGFSRACYIDTNVLH]A

[SPAT][R]



TCCG

GC content: 75 %.

Can form the following amino acids in 6 translation frames:

R[RSIKTNM]

[PGFSRACYIDTNVLH]G

[SPAT][ED]

S[GEADV]

[PGFSRACYIDTNVLH]P

[IVFL][R]



CGGG

GC content: 100 %.

Has 3 identical bases in a row. However, this has not shown to be very important.

Can form the following amino acids in 6 translation frames:

P[GEADV]
[PGFSRACYIDTNVLH]P
[SPAT][R]
R[GEADV]
[PGFSRACYIDTNVLH]G
[SPAT][G]



TACG

GC content: 50 %.

Can form the following amino acids in 6 translation frames:

R[RSIKTNM]

[PGFSRACYIDTNVLH]V

[SPAT][Y*]

Y[GEADV]

[PGFSRACYIDTNVLH]T

[IVL][R]



CGTT

GC content: 50 %.

Can form the following amino acids in 6 translation frames:

N[GEADV]
[PEGRSA*IKTQVL]T
[EQK*][R]
R[FWSCY*L]
[PGFSRACYIDTNVLH]V
[SPAT][LF]



GTAG

GC content: 50 %.

Can form the following amino acids in 6 translation frames:

L[PRQLH]
[PGFSRACYIDTNVLH]Y
[SPAT][T]
V[GEADV]
[PEGWRSA*KTQMVL]*
[SGCR][RS]

Page 63 of 138



CTAG

GC content: 50 %.

The overhang is palindromic, cannot be used for DNA assembly.

Can form the following amino acids in 6 translation frames:

L[GEADV]

[PGFSRACYIDTNVLH]*

[SPAT][RS]

L[GEADV]

[PGFSRACYIDTNVLH]*

[SPAT][RS]



CTAT

GC content: 25 %.

[SPAT][MI]

Can form the following amino acids in 6 translation frames:

I[GEADV]
[PEGRSA*IKTQVL]*
[HDYN][RS]
L[FWSCY*L]
[PGFSRACYIDTNVLH]Y

Page 65 of 138



GGAG

GC content: 75 %.

Can form the following amino acids in 6 translation frames:

L[PRQLH]

[PGFSRACYIDTNVLH]S

[SPAT][P]

G[GEADV]

[PEGWRSA*KTQMVL]E

[RGW][RS]



CTCG

GC content: 75 %.

Can form the following amino acids in 6 translation frames:

R[GEADV]

[PGFSRACYIDTNVLH]E

[SPAT][RS]

L[GEADV]

[PGFSRACYIDTNVLH]S

[SPAT][R]



CTCT

GC content: 50 %.

[SPAT][L]

Can form the following amino acids in 6 translation frames:

R[GEADV]
[PEGRSA*IKTQVL]E
[EQK*][RS]
L[FWSCY*L]
[PGFSRACYIDTNVLH]S

Page 68 of 138



CTGG

GC content: 75 %.

Can form the following amino acids in 6 translation frames:

P[GEADV]

[PGFSRACYIDTNVLH]Q

[SPAT][RS]

L[GEADV]

[PGFSRACYIDTNVLH]W

[SPAT][G]



GAAG

GC content: 50 %.

Can form the following amino acids in 6 translation frames:

L[PRQLH]

[PGFSRACYIDTNVLH]F

[SPAT][S]

E[GEADV]

[PEGWRSA*KTQMVL]K

[RG*][RS]



CTTG

GC content: 50 %.

Can form the following amino acids in 6 translation frames:

Q[GEADV]

[PGFSRACYIDTNVLH]K

[SPAT][RS]

L[GEADV]

[PGFSRACYIDTNVLH]L

[SPAT][WC*]





GC content: 25 %.

Has 3 identical bases in a row. However, this has not shown to be very important.

Can form the following amino acids in 6 translation frames:

K[GEADV]
[PEGRSA*IKTQVL]K
[EQK*][RS]
L[FWSCY*L]
[PGFSRACYIDTNVLH]F
[SPAT][LF]



TTTC

GC content: 25 %.

Has 3 identical bases in a row. However, this has not shown to be very important.

Can form the following amino acids in 6 translation frames:

E[RSIKTNM]
[PEGWRSA*KTQMVL]K
[RG*][KN]
F[PRQLH]
[PGFSRACYIDTNVLH]F
[IVFL][S]



GTTC

GC content: 50 %.

Can form the following amino acids in 6 translation frames:

E[PRQLH]

[PEGWRSA*KTQMVL]N

[RG*][T]

V[PRQLH]

[PEGWRSA*KTQMVL]F

[SGCR][S]



GGTC

GC content: 75 %.

[RGW][S]

Can form the following amino acids in 6 translation frames:

D[PRQLH]
[PEGWRSA*KTQMVL]T
[RG*][P]
G[PRQLH]
[PEGWRSA*KTQMVL]V



GACG

GC content: 75 %.

Can form the following amino acids in 6 translation frames:

R[PRQLH]

[PGFSRACYIDTNVLH]V

[SPAT][S]

D[GEADV]

[PEGWRSA*KTQMVL]T

[RG*][R]



GACT

GC content: 50 %.

Can form the following amino acids in 6 translation frames:

S[PRQLH]
[PEGRSA*IKTQVL]V
[EQK*][S]
D[FWSCY*L]
[PEGWRSA*KTQMVL]T
[RG*][L]



GCTC

GC content: 75 %.

Can form the following amino acids in 6 translation frames:

E[PRQLH]
[PEGWRSA*KTQMVL]S
[RG*][A]
A[PRQLH]

[PEGWRSA*KTQMVL]L

[SGCR][S]



GAGT

GC content: 50 %.

Can form the following amino acids in 6 translation frames:

T[PRQLH]

[PEGRSA*IKTQVL]L

[HDYN][S]

E[FWSCY*L]

[PEGWRSA*KTQMVL]S

[RG*][V]



GATC

GC content: 50 %.

The overhang is palindromic, cannot be used for DNA assembly.

Can form the following amino acids in 6 translation frames:

D[PRQLH]

[PEGWRSA*KTQMVL]I

[RG*][S]

D[PRQLH]

[PEGWRSA*KTQMVL]I

[RG*][S]



TTGC

GC content: 50 %.

[IVFL][A]

Can form the following amino acids in 6 translation frames:

A[RSIKTNM]
[PEGWRSA*KTQMVL]Q
[SGCR][KN]
L[PRQLH]
[PGFSRACYIDTNVLH]C

Page 81 of 138



GCAG

GC content: 75 %.

Can form the following amino acids in 6 translation frames:

L[PRQLH] [PGFSRACYIDTNVLH]C

[SPAT][A] A[GEADV]

[PEGWRSA*KTQMVL]Q

[SGCR][RS]



GCAT

GC content: 50 %.

The overhang contains the start codon ATG.

Can form the following amino acids in 6 translation frames:

M[PRQLH]
[PEGRSA*IKTQVL]C
[HDYN][A]
A[FWSCY*L]

[PEGWRSA*KTQMVL]H

[SGCR][MI]



GGGC

GC content: 100 %.

Has 3 identical bases in a row. However, this has not shown to be very important.

Can form the following amino acids in 6 translation frames:

A[PRQLH]
[PEGWRSA*KTQMVL]P
[SGCR][P]
G[PRQLH]
[PEGWRSA*KTQMVL]G

[RGW][A]



✓ CGGC

GCCG

GC content: 100 %.

Can form the following amino acids in 6 translation frames:

R[PRQLH]

[PGFSRACYIDTNVLH]G

[SPAT][A]

A[GEADV]

[PEGWRSA*KTQMVL]P

[SGCR][R]



GCCT

GC content: 75 %.

Can form the following amino acids in 6 translation frames:

R[PRQLH] [PEGRSA*IKTQVL]G [EQK*][A] A[FWSCY*L] [PEGWRSA*KTQMVL]P [SGCR][L]



TCGC

GC content: 75 %.

Can form the following amino acids in 6 translation frames:

A[RSIKTNM]

[PEGWRSA*KTQMVL]R

[SGCR][ED]

S[PRQLH]

[PGFSRACYIDTNVLH]R

[IVFL][A]



X GCGC

GCGC

GC content: 100 %.

The overhang is palindromic, cannot be used for DNA assembly.

Can form the following amino acids in 6 translation frames:

A[PRQLH]

[PEGWRSA*KTQMVL]R

[SGCR][A]

A[PRQLH]

[PEGWRSA*KTQMVL]R

[SGCR][A]



TAGC

GC content: 50 %.

[IVL][A]

Can form the following amino acids in 6 translation frames:

A[RSIKTNM]
[PEGWRSA*KTQMVL]L
[SGCR][Y*]
*[PRQLH]
[PGFSRACYIDTNVLH]S



GCTG

GC content: 75 %.

Can form the following amino acids in 6 translation frames:

Q[PRQLH]
[PGFSRACYIDTNVLH]S
[SPAT][A]
A[GEADV]
[PEGWRSA*KTQMVL]L
[SGCR][WC*]



GCTT

GC content: 50 %.

Can form the following amino acids in 6 translation frames:

K[PRQLH]
[PEGRSA*IKTQVL]S
[EQK*][A]
A[FWSCY*L]
[PEGWRSA*KTQMVL]L
[SGCR][LF]



GTCC

GC content: 75 %.

Can form the following amino acids in 6 translation frames:

G[PRQLH]

[PEGWRSA*KTQMVL]D

[RGW][T]

V[PRQLH]

[PEGWRSA*KTQMVL]S

[SGCR][P]



X GGCC

GGCC

GC content: 100 %.

The overhang is palindromic, cannot be used for DNA assembly.

Can form the following amino acids in 6 translation frames:

G[PRQLH]

[PEGWRSA*KTQMVL]A

[RGW][P]

G[PRQLH]

[PEGWRSA*KTQMVL]A

[RGW][P]





GC content: 75 %.

Has 3 identical bases in a row. However, this has not shown to be very important.

Can form the following amino acids in 6 translation frames:

G[RSIKTNM]
[PEGWRSA*KTQMVL]G
[RGW][ED]
S[PRQLH]
[PGFSRACYIDTNVLH]P
[IVFL][P]



GGGG

GC content: 100 %.

Has 3 identical bases in a row. However, this has not shown to be very important.

Can form the following amino acids in 6 translation frames:

P[PRQLH]
[PGFSRACYIDTNVLH]P
[SPAT][P]
G[GEADV]
[PEGWRSA*KTQMVL]G
[RGW][G]



GGGT

GC content: 75 %.

Has 3 identical bases in a row. However, this has not shown to be very important.

Can form the following amino acids in 6 translation frames:

T[PRQLH]
[PEGRSA*IKTQVL]P
[HDYN][P]
G[FWSCY*L]
[PEGWRSA*KTQMVL]G
[RGW][V]



TTAC

GC content: 25 %.

Can form the following amino acids in 6 translation frames:

V[RSIKTNM]

[PEGWRSA*KTQMVL]*

[SGCR][KN]

L[PRQLH]

[PGFSRACYIDTNVLH]Y

[IVFL][T]



GTAC

GC content: 50 %.

The overhang is palindromic, cannot be used for DNA assembly.

Can form the following amino acids in 6 translation frames:

V[PRQLH]

[PEGWRSA*KTQMVL]Y

[SGCR][T]

V[PRQLH]

[PEGWRSA*KTQMVL]Y

[SGCR][T]



GTAT

GC content: 25 %.

Can form the following amino acids in 6 translation frames:

I[PRQLH]
[PEGRSA*IKTQVL]Y
[HDYN][T]
V[FWSCY*L]
[PEGWRSA*KTQMVL]Y
[SGCR][MI]



TGAC

GC content: 50 %.

The overhang contains a stop codon (TAA, TAG or TGA).

Can form the following amino acids in 6 translation frames:

V[RSIKTNM]
[PEGWRSA*KTQMVL]S
[SGCR][HQ]
*[PRQLH]
[PGFSRACYIDTNVLH]D
[MVL][T]



GTCT

GC content: 50 %.

Can form the following amino acids in 6 translation frames:

R[PRQLH]
[PEGRSA*IKTQVL]D
[EQK*][T]
V[FWSCY*L]
[PEGWRSA*KTQMVL]S
[SGCR][L]



TCAC

GC content: 50 %.

[IVFL][T]

The overhang contains a stop codon (TAA, TAG or TGA).

Can form the following amino acids in 6 translation frames:

V[RSIKTNM]
[PEGWRSA*KTQMVL]*
[SGCR][ED]
S[PRQLH]
[PGFSRACYIDTNVLH]H



GTGC

GC content: 75 %.

Can form the following amino acids in 6 translation frames:

A[PRQLH]

[PEGWRSA*KTQMVL]H

[SGCR][T]

V[PRQLH]

[PEGWRSA*KTQMVL]C

[SGCR][A]



GTGT

GC content: 50 %.

Can form the following amino acids in 6 translation frames:

T[PRQLH]
[PEGRSA*IKTQVL]H
[HDYN][T]
V[FWSCY*L]
[PEGWRSA*KTQMVL]C

[SGCR][V]



TAAC

GC content: 25 %.

Can form the following amino acids in 6 translation frames:

V[RSIKTNM]

[PEGWRSA*KTQMVL]L

[SGCR][Y*]

*[PRQLH]

[PGFSRACYIDTNVLH]N

[IVL][T]



TAAG

GC content: 25 %.

Can form the following amino acids in 6 translation frames:

L[RSIKTNM]

[PGFSRACYIDTNVLH]L

[SPAT][Y*]

*[GEADV]

[PGFSRACYIDTNVLH]K

[IVL][RS]



TACC

GC content: 50 %.

Can form the following amino acids in 6 translation frames:

G[RSIKTNM]

[PEGWRSA*KTQMVL]V

 $[RGW][Y^*]$

Y[PRQLH]

[PGFSRACYIDTNVLH]T

[IVL][P]



TAGG

GC content: 50 %.

Can form the following amino acids in 6 translation frames:

P[RSIKTNM]

[PGFSRACYIDTNVLH]L

[SPAT][Y*]

*[GEADV]

[PGFSRACYIDTNVLH]R

[IVL][G]



TATA

GC content: 0 %.

The overhang is palindromic, cannot be used for DNA assembly.

Can form the following amino acids in 6 translation frames:

Y[RSIKTNM]
[PGFSRACYIDTNVLH]I
[IVL][Y*]
Y[RSIKTNM]
[PGFSRACYIDTNVLH]I
[IVL][Y*]



TATC

GC content: 25 %.

Can form the following amino acids in 6 translation frames:

D[RSIKTNM]

[PEGWRSA*KTQMVL]I

[RG*][Y*]

Y[PRQLH]

[PGFSRACYIDTNVLH]I

[IVL][S]



TATG

GC content: 25 %.

The overhang contains the start codon ATG.

Can form the following amino acids in 6 translation frames:

H[RSIKTNM]

[PGFSRACYIDTNVLH]I

[SPAT][Y*]

Y[GEADV]

[PGFSRACYIDTNVLH]M

[IVL][WC*]



TTGA

GC content: 25 %.

The overhang contains a stop codon (TAA, TAG or TGA).

Can form the following amino acids in 6 translation frames:

S[RSIKTNM]

[PGFSRACYIDTNVLH]Q

[IVFL][KN]

L[RSIKTNM]

[PGFSRACYIDTNVLH]*

[IVFL][ED]



TCAG

GC content: 50 %.

The overhang contains a stop codon (TAA, TAG or TGA).

Can form the following amino acids in 6 translation frames:

L[RSIKTNM]
[PGFSRACYIDTNVLH]*
[SPAT][ED]
S[GEADV]
[PGFSRACYIDTNVLH]Q
[IVFL][RS]



TCGA

GC content: 50 %.

The overhang is palindromic, cannot be used for DNA assembly.

Can form the following amino acids in 6 translation frames:

S[RSIKTNM]

[PGFSRACYIDTNVLH]R

[IVFL][ED]

S[RSIKTNM]

[PGFSRACYIDTNVLH]R

[IVFL][ED]



TCGG

GC content: 75 %.

Can form the following amino acids in 6 translation frames:

P[RSIKTNM]

[PGFSRACYIDTNVLH]R

[SPAT][ED]

S[GEADV]

[PGFSRACYIDTNVLH]R

[IVFL][G]



TCTA

GC content: 25 %.

[IVFL][Y*]

Can form the following amino acids in 6 translation frames:

*[RSIKTNM]
[PGFSRACYIDTNVLH]R
[IVL][ED]
S[RSIKTNM]
[PGFSRACYIDTNVLH]L



TCTC

GC content: 50 %.

Can form the following amino acids in 6 translation frames:

E[RSIKTNM] [PEGWRSA*KTQMVL]R [RG*][ED] S[PRQLH]

[PGFSRACYIDTNVLH]L

[IVFL][S]



TCTT

GC content: 25 %.

Can form the following amino acids in 6 translation frames:

K[RSIKTNM]

[PEGRSA*IKTQVL]R

[EQK*][ED]

S[FWSCY*L]

[PGFSRACYIDTNVLH]L

[IVFL][LF]



TTCA

GC content: 25 %.

The overhang contains a stop codon (TAA, TAG or TGA).

Can form the following amino acids in 6 translation frames:

*[RSIKTNM]

[PGFSRACYIDTNVLH]E

[MVL][KN]

F[RSIKTNM]

[PGFSRACYIDTNVLH]S

[IVFL][HQ]



TGAG

GC content: 50 %.

The overhang contains a stop codon (TAA, TAG or TGA).

Can form the following amino acids in 6 translation frames:

L[RSIKTNM]
[PGFSRACYIDTNVLH]S
[SPAT][HQ]
*[GEADV]
[PGFSRACYIDTNVLH]E
[MVL][RS]



TGAT

GC content: 25 %.

The overhang contains a stop codon (TAA, TAG or TGA).

Can form the following amino acids in 6 translation frames:

I[RSIKTNM] [PEGRSA*IKTQVL]S

[HDYN][HQ]

*[FWSCY*L]

[PGFSRACYIDTNVLH]D

[MVL][MI]



TGCA

GC content: 50 %.

The overhang is palindromic, cannot be used for DNA assembly.

Can form the following amino acids in 6 translation frames:

C[RSIKTNM]
[PGFSRACYIDTNVLH]A
[MVL][HQ]
C[RSIKTNM]
[PGFSRACYIDTNVLH]A
[MVL][HQ]



TGCC

GC content: 75 %.

Can form the following amino acids in 6 translation frames:

G[RSIKTNM]
[PEGWRSA*KTQMVL]A
[RGW][HQ]
C[PRQLH]

[PGFSRACYIDTNVLH]A

[MVL][P]



TGCT

GC content: 50 %.

Can form the following amino acids in 6 translation frames:

S[RSIKTNM]
[PEGRSA*IKTQVL]A
[EQK*][HQ]
C[FWSCY*L]
[PGFSRACYIDTNVLH]A

[MVL][L]



TGGA

GC content: 50 %.

[MVL][ED]

Can form the following amino acids in 6 translation frames:

S[RSIKTNM]
[PGFSRACYIDTNVLH]P
[IVFL][HQ]
W[RSIKTNM]
[PGFSRACYIDTNVLH]G



TGGC

GC content: 75 %.

Can form the following amino acids in 6 translation frames:

A[RSIKTNM]
[PEGWRSA*KTQMVL]P
[SGCR][HQ]
W[PRQLH]
[PGFSRACYIDTNVLH]G
[MVL][A]



TGTA

GC content: 25 %.

Can form the following amino acids in 6 translation frames:

Y[RSIKTNM]
[PGFSRACYIDTNVLH]T
[IVL][HQ]
C[RSIKTNM]
[PGFSRACYIDTNVLH]V
[MVL][Y*]



TGTC

GC content: 50 %.

Can form the following amino acids in 6 translation frames:

D[RSIKTNM]

[PEGWRSA*KTQMVL]T

[RG*][HQ]

C[PRQLH]

[PGFSRACYIDTNVLH]V

[MVL][S]





GC content: 0 %.

The overhang is palindromic, cannot be used for DNA assembly.

Can form the following amino acids in 6 translation frames:

L[RSIKTNM]
[PGFSRACYIDTNVLH]*
[IVFL][KN]
L[RSIKTNM]
[PGFSRACYIDTNVLH]*
[IVFL][KN]



TTAG

GC content: 25 %.

Can form the following amino acids in 6 translation frames:

L[RSIKTNM]

[PGFSRACYIDTNVLH]*

[SPAT][KN]

L[GEADV]

[PGFSRACYIDTNVLH]*

[IVFL][RS]



TTAT

GC content: 0 %.

Can form the following amino acids in 6 translation frames:

I[RSIKTNM]

[PEGRSA*IKTQVL]*

[HDYN][KN]

L[FWSCY*L]

[PGFSRACYIDTNVLH]Y

[IVFL][MI]



TTCC

GC content: 50 %.

Can form the following amino acids in 6 translation frames:

G[RSIKTNM]

[PEGWRSA*KTQMVL]E

[RGW][KN]

F[PRQLH]

[PGFSRACYIDTNVLH]S

[IVFL][P]



TTCG

GC content: 50 %.

Can form the following amino acids in 6 translation frames:

R[RSIKTNM]

[PGFSRACYIDTNVLH]E

[SPAT][KN]

F[GEADV]

[PGFSRACYIDTNVLH]S

[IVFL][R]



TTCT

GC content: 25 %.

Can form the following amino acids in 6 translation frames:

R[RSIKTNM]

[PEGRSA*IKTQVL]E

[EQK*][KN]

F[FWSCY*L]

[PGFSRACYIDTNVLH]S

[IVFL][L]



TTGG

GC content: 50 %.

Can form the following amino acids in 6 translation frames:

P[RSIKTNM]

[PGFSRACYIDTNVLH]Q

[SPAT][KN]

L[GEADV]

[PGFSRACYIDTNVLH]W

[IVFL][G]





GC content: 0 %.

Has 3 identical bases in a row. However, this has not shown to be very important.

Can form the following amino acids in 6 translation frames:

*[RSIKTNM]
[PGFSRACYIDTNVLH]K
[IVL][KN]
F[RSIKTNM]
[PGFSRACYIDTNVLH]L
[IVFL][Y*]





GC content: 0 %.

Has 3 identical bases in a row. However, this has not shown to be very important.

Can form the following amino acids in 6 translation frames:

K[RSIKTNM]
[PEGRSA*IKTQVL]K
[EQK*][KN]
F[FWSCY*L]
[PGFSRACYIDTNVLH]F
[IVFL][LF]

Appendix

The report consists of 3 sections: summary, overhangs, appendix.

Summary page(s)

The first page summarises the compendium.

Overhang pages

Each overhang is analysed separately. The result is summarised with a symbol:

☑: good overhang

Overhangs are unpaired nucleotides at the end of a double-stranded linear DNA molecule. Overhangs can be on either strand; 5' or 3' overhangs. After DNA ligation with another DNA with a complementary overhang, these remain in the sequence as fusion sites ("scars").

Overhang sets

Use the GoldenHinges Python package to generate a set of mutually compatible overhangs that can be used for DNA assembly.